



Recombinant Human Olfactory receptor 2J3 (OR2J3)

Product Code	CSB-CF016583HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O76001
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNDDGKVNASSEGYPILVGFSNWPHELVVIFVVVLIFYLMTLIGNLFIILSYLDS HLHT PMYFFLSNLSFLDLCYTTSSIPQLLVNLWGPEKTISYAGCMIQLYFVLALGTTEC VLLVV MSYDRYAAVCRPLHYTVLMHPRFCHLLAVASWVSGFTNSALHSSFTFWVPLC GHRQVDHF FCEVPALLRLSCVDTHVNELTLMITSSIFVLIPLILILTSYGAIVRAVLRMQSTTGL QKV FGTCGAHLMAVSLFFIPAMCIYLQPPSGNSQDQGKFIALFYTVVTPSLNPLIYTL RNKVV RGAVKRLMGWE
Source	in vitro E.coli expression system
Target Names	OR2J3
Protein Names	Recommended name: Olfactory receptor 2J3 Alternative name(s): Hs6M1-3 Olfactory receptor OR6-16 Short name= OR6-6 Short name= Olfactory receptor 6-6
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.