



Recombinant Human Olfactory receptor 6B1 (OR6B1)

Product Code	CSB-CF016994HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95007
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MELENQTRVTKFILVGFPGLSMRAAMFLIFLVAYILTVAENVIIILLVLQNRPLHK PMY FFLANLSFLETWYISVTVPKLLFSFWSVNNSISFTLCMIQLYFFIALMCTECVLLA AMAY DRYVAICRPLHYPTIMSHGLCFRLALGSWAIGFGISLAKIYFISCLSFCGPNVINH FFCD ISPVLNLSCTDMSITELVDFILALVIFLPLFITVLSYGCILATILCMPTGKQKAFST CA SHLVVVTIFYSAIIFMYARPRVIHAFNMNKIISIFYAIVTPSLNPFYICLRNREVKEA LK KLAYCQASRSD
Source	in vitro E.coli expression system
Target Names	OR6B1
Protein Names	Recommended name: Olfactory receptor 6B1 Alternative name(s): Olfactory receptor 7-3 Short name= OR7-3 Olfactory receptor OR7-9
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.