



Recombinant Human Olfactory receptor 4F6 (OR4F6)

Product Code	CSB-CF850886HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGB9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDEANHSVSEFVFLGLSDSRKIQLLLFLFFSVFYVSSLMGNLLIVLTVTSDPRL QSPMY FLLANLSIINLVFCSSTAPKMIYDLFRKHKTISFGGCVVQIFFIHAVGGTEMVLLIA MAF DRYVAICKPLHYLTIMNPQRCILFLVISWIIGIIHSVIQLAFVVDLLFCGPNELDSFF CD LPRFIKLACIETYTLGFMVTANSGFISLASFLILIISYIFILVTVQKKSSGGIFKAFSM L SAHVIVVVLVFGPLIFFYIFPFPTSHLDKFLAIFDAVITPVLNPNVIYTFRNKEMMVA MRR RCSQFVNYSKIF
Source	in vitro E.coli expression system
Target Names	OR4F6
Protein Names	Recommended name: Olfactory receptor 4F6 Alternative name(s): Olfactory receptor 4F12 Olfactory receptor OR15-15
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.