



Recombinant Human Olfactory receptor 8K3 (OR8K3)

Product Code	CSB-CF854151HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH51
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEQHNLTTVNEFILTGITDIAELQAPLFALFLMIYVISVMGNLGMIVLTKLDSRLQ TPMY FFLRHLAFMDLGYSTTVGPKMLVNFVVDKNIISYYFCATQLAFFLVFIGSELFILS AMSY DLYVAICNPLLYTVIMSRRVCQVLVAIPYLYCTFISLLVTIKIFTLSFCGYNVISHFY CD SLPLLPLLCNTHIEIILIFAIDLISLLIVLLSYLLILVAILRMNSAGRQKAFSTC GAHLTVVIVFYGTLLFMVYVQPKSSHSDTKVASIFYTLVIPMLNPLIYSLRNKD VKYAL RRTWNNLCNIFV
Source	in vitro E.coli expression system
Target Names	OR8K3
Protein Names	Recommended name: Olfactory receptor 8K3 Alternative name(s): Olfactory receptor OR11-181
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.