



Recombinant Human Olfactory receptor 7D2 (OR7D2)

Product Code	CSB-CF857016HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96RA2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEAGNQTGFLEFILLGLSEDPELQPFIFGLFLSMYLVTVLGNLLIILAISSDSHLH TPMY FFLSNLSWVDICFSTCIVPKMLVNIQTENKAISYMDCLTQVYFSMFFPILDTLLLT VMAY DRFVAVCHPLHYMIIMNPHLCGLLVFVTWLIGVMTSLLHISLMMHLIFCKDFEIP HFFCE LTYILQLACSDTFLNSTLIYFMTGVLGVFPLLGIIFSYSRIASSIRKMSSSGGKQK ALST CGSHLSVVSIFYGTGIGVHFTSAVTHSSQKISVASVMYTVVTPMLNPFYSLRN KDVKGA LGSLLSRAASCL
Source	in vitro E.coli expression system
Target Names	OR7D2
Protein Names	Recommended name: Olfactory receptor 7D2 Alternative name(s): HTPCRH03 Olfactory receptor 19-4 Short name= OR19-4 Olfactory receptor OR19-10
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.