



Recombinant Human Olfactory receptor 51B6 (OR51B6)

Product Code	CSB-CF875662HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H340
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGLNKSASTFQLTGFPMEKAHHWIFIPLLAAYISILLGNGTLLFLIRNDHNLHE PMYYF LAMLAATDLGVTLTTMPTVLGVWLWDHREIGHGACFSQAYFIHTLSVMESGVLL AMAYDC FITIRSPLRYTSILTNTQVMKIGVRVLTRAGLSIMPIVRLHWFPYCRSHVLSHAF CLHQ DVIKLACADITFNRLYPVVVLFAMVLLDFLIIFFSYILILKTVMGIGSGGERAKALN TCV SHICCILVFYVTVVCLTFIHRFGKHVPHVVHITMSYIHFLFPPFMNPFYISIKTKQI QSG ILRLFSLPHSRA
Source	in vitro E.coli expression system
Target Names	OR51B6
Protein Names	Recommended name: Olfactory receptor 51B6 Alternative name(s): Odorant receptor HOR5'beta6
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.