



Recombinant Human Olfactory receptor 2J2 (OR2J2)

Product Code	CSB-CF016582HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O76002
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMIKKNASSEDFILLGFSNWPQLEVVLFVVILIFYLMTLTGNLFIILSYVDSHLH TPM YFFLSNLSFLDLCHTTSSIPQLLVNLRGPEKTISYAGCMVQLYFVLALGIAECVL LVVMS YDRYVAVCRPLHYTVLMHPRFCHLLAAASWVIGFTISALHSSFTFWVPLCGHR LVDHFFC EVPALLRLSCVDTHANELTLMVMSSIFVLIPLILILTAYGAIARAVLSMQSTTGLQ KVFR TCGAHLMVVSLFFIPVMCMYLQPPSENSPDQGKFIALFYTVVTPSLNPLIYTLR NKHVKG AAKRLLGWEWGK
Source	in vitro E.coli expression system
Target Names	OR2J2
Protein Names	Recommended name: Olfactory receptor 2J2 Alternative name(s): Hs6M1-6 Olfactory receptor 6-8 Short name= OR6-8 Olfactory receptor OR6-19
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.