



Recombinant Human Olfactory receptor 6V1 (OR6V1)

Product Code	CSB-CF839777HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N148
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MANLSQPSEFVLLGFSSFGELQALLYGPFLMLYLLAFMGNTIIIVMVIADTHLHT PMYFF LGNFSLLEILVTMTAVPRMLSDLLVPHKVITFTGCMVQFYFHFSLGSTSFLILTD MALDR FVAICHPLRYGTLMSCRAMCVQLAGAAWAAPFLAMVPTVLSRAHLDYCHGDVIN HFFCDNE PLLQLSCSDTRLLEFWDFLMALTFVLSSFLVTLISYGYIVTTVLRIPSASSCQKAF STCG SHLTLVFIGYSSTIFLYVRPGKAHSVQVRKVVALVTSVLTPLNPFILTFCNQTV KTVLQ GQMQRLLKGLCKAQ
Source	in vitro E.coli expression system
Target Names	OR6V1
Protein Names	Recommended name: Olfactory receptor 6V1 Alternative name(s): Olfactory receptor OR7-3
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.