



Recombinant Human Olfactory receptor 4Q3 (OR4Q3)

Product Code	CSB-CF839907HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH05
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKKEQDSNVTEFVLLGLSSSWELQLFLFLLFLFFYIAIVLGNLLIVVTVQAAHLL QSPM YYFLGHLSFIDLCLSCVTPKMLGDFLQQGKSISFSGCLAQIYFLHFLGASEMFL LTVMA YDRYVAICNPLRYLTVMNPQLCLWLVLACWCGGFIHSIMQVILVIQLPFCGPNE LDFYC DVPQVIKLACMDTYVVEVLVIANSGLLSLVCFLVLLFSYAILITLRTHFCQGQNK VFST CASHLTVVSLIFVPCVFIYLRPFCSFSVDKIFSLFYTVITPMLNPLIYTLRNTDMKT AMK KLRIKPCGIPLPC
Source	in vitro E.coli expression system
Target Names	OR4Q3
Protein Names	Recommended name: Olfactory receptor 4Q3 Alternative name(s): Olfactory receptor 4Q4 Olfactory receptor OR14-3
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.