



Recombinant Human Olfactory receptor 8B3 (OR8B3)

Product Code	CSB-CF844061HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGG8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLARNNSLVTEFILAGLTDHPEFQQPLFFLFLVVYIVTMVGNLGLIILFGLNSHLH TPMY YFLFNLSFIDLCYSSVFTPKMLMNFVSKKNIISYVGCMTQLFFFLFFVISECYMLT SMAY DRYVAICNPLLYKVTMSHQVCSMLTFAAYIMGLAGATAHTGCMLRLTFCSANII NHYLCD ILPLLQLSCTSTYVNEVVVLIVVGINIMVPSCILISYVFIVTSILHIKSTQGRSKAF ST CSSHVIALSLFFGSAAFMYIKYSSGSMEQGKVSSVFYTNVPMNLPLIYSLRNK DVKVAL RKALIKIQRRNIF
Source	in vitro E.coli expression system
Target Names	OR8B3
Protein Names	Recommended name: Olfactory receptor 8B3 Alternative name(s): Olfactory receptor OR11-311
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.