



Recombinant Human Olfactory receptor 4M2 (OR4M2)

Product Code	CSB-CF847716HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGB6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METANYTKVTEFVLTGLSQTPEVQLVLFVIFLSFYLFILPGNIIICTISLDPHLTSP MY FLLANLAFLDIWYSSITAPEMLIDFFVERKIIISFDGCIQLFFLHFAGASEMFLTV MAF DLYTAICRPLHYATIMNQRLCCILVALSWRGGFIHSIIQVALIVRLPFCGPNELDS YFCD ITQVVRIACANTFPEELVMICSSGLISVVCLIALLLMSYAFLLALFKKLSGSGENTN RAMS TCYSHITIVLMFGPSIYIYARPFDSFSLDKVVSFNTLIFPLRNPIIYTLRNKEVK AAM RKLVTKYILCKEK
Source	in vitro E.coli expression system
Target Names	OR4M2
Protein Names	Recommended name: Olfactory receptor 4M2 Alternative name(s): Olfactory receptor OR15-3
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.