



Recombinant Human Olfactory receptor 51A4 (OR51A4)

Product Code	CSB-CF847724HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGJ6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSIINTSYVEITTFVLVGMPGLEIAHIWISIPICSMYLIAILGNGTILFIIKTEPSLHE P MYYFLSMLAMSDLGLSLSSLPTVLSIFLFNAPEISSNACFAQEFFIHGFSVLESS VLLIM SFDRFLAIHNPLRYTSILTTRVAQIGIVFSFKSMMLVLPFPFTLRNLRVCKKNQL SHSY CLHQDVMKLACSDNRIDVIYGFFGALCLMVDFILIAVSYTLILKTVLGIASKKEQL KALN TCVSHICAVIIFYLPILAVVHRFARHVSPLINVLMANVLLLVPPLTNPIVYCVKT KQI RVRVVAKLCQRKI
Source	in vitro E.coli expression system
Target Names	OR51A4
Protein Names	Recommended name: Olfactory receptor 51A4
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.