



Recombinant Human Olfactory receptor 2B6 (OR2B6)

Product Code	CSB-CF016561HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58173
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNWWNDSIIQEFILLGFSDRPWLEFLLVVFLISYTVTIFGNLTIIILVSRLDTKLHT PMY FFLTNLSLLDLCYTTCTVPQMLVNLCSIRKVISYRGCVAQLFIFLALGATEYLLLA VMSF DRFVAICRPLHYSVIMHQRLCLQLAAASWVTGFSNSVWLSTLTQLPLCDPYVI DHFLCE VPALLKLSCVETTANEAEFLVSELFHLIPLTLILISYAFIVRAVLRIQSAEGRQKA FGT CGSHLIVVSLFYSTAVSVYLVQPPSPSSKDQGKMVSLFYGIIAPMLNPLIYTLRNK EVKEG FKRLVARVFLIKK
Source	in vitro E.coli expression system
Target Names	OR2B6
Protein Names	Recommended name: Olfactory receptor 2B6 Alternative name(s): Hs6M1-32 Olfactory receptor 2B1 Olfactory receptor 2B5 Olfactory receptor 5-40 Short name= OR5-40 Olfactory receptor 6-31 Short name= OR6-31 Olfactory
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.