



Recombinant Human Olfactory receptor 5H1 (OR5H1)

Product Code	CSB-CF016946HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NKK0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEEENATLLTEFVLTGFLYQPQWKIPLFLAFLVIYLITIMGNLGLIAVIWKDPHLHI PMY LLLGNLAFVDAWISSTVTPKMLNNFLAKSKMISLSECKIQFFSFAISVTTECFLLA TMAY DRYVAICKPLLYPAIMTNGLCIRLLILSYVGGILHALIHEGFLFRLTFCNSNIVHHIY CD TIPLSKISCTDSSINFLMVFIFSGSIQVFSIVTILVSYTFVLFVFAILKKKSDKGVRKAF ST CGAHLFSVSLYYGPLLFIYVGPASPQADDQDMVEPLFYTVIIPLLNPIIYSLRNKQ VTVS FTKMLKKHVKVS
Source	in vitro E.coli expression system
Target Names	OR5H1
Protein Names	Recommended name: Olfactory receptor 5H1 Alternative name(s): HTPCRX14
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.