



Recombinant Human Olfactory receptor 4D6 (OR4D6)

Product Code	CSB-CF818797HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGJ1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDQINHTNVKEFFFLLELTRSRELEFFLVVFFAVYVATVLGNALIVVTITCESRL HTPMY FLLRNKSVLDIVFSSITVPKFLVDLLSDRKTISYNDCMAQIFFFHFAGGADIFFLS VMAY DRYLAIKPLHYVTMMRKEVWVALVVASWVSGGLHSIIQVILMLPFPFCGPNTL DAFYCY VLQVVKLACTDTFALELFMISNGLVTLWFLLLLGSYTVILVMLRSHSGEGRN KALSTC TSHMLVVTLHFVPCVYIYCRPFMTLPMDTTISINNTVITPMLNPIIYSLRNQEMKS AMQR LQRRLGPSESARKWG
Source	in vitro E.coli expression system
Target Names	OR4D6
Protein Names	Recommended name: Olfactory receptor 4D6 Alternative name(s): Olfactory receptor OR11-250
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.