



Recombinant Human Olfactory receptor 4K2 (OR4K2)

Product Code	CSB-CF850887HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGD2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDVGNKSTMSEFVLLGLSNSWELQMFFFMVFSLLYVATMVGNSLIVITVIVDPH LHSPMY FLLTNLSIIDMSLASFATPKMITDYLTGHKTISFDGCLTQIFFLHLFTGTEIILLMAM SF DRYIAICKPLHYASVISPVQVCVALVVASWIMGVMHMSMSQVIFALTLPCGPEYEV DSFFCD LPVVFQLACVDITYVLGLFMISTSGIIALSCFIVLFNSYVIVLVTVKHHSSRGSSKA LSTC TAHFIVVFLFFGPCIFIYMWPLSSFLTDKILSVFYTIFTPTLNPIIYTLRNQEVKIAM RK LKNRFLNFKAMPS
Source	in vitro E.coli expression system
Target Names	OR4K2
Protein Names	Recommended name: Olfactory receptor 4K2 Alternative name(s): Olfactory receptor OR14-15
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.