



Recombinant Human Olfactory receptor 2T5 (OR2T5)

Product Code	CSB-CF740230HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6IEZ7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MANITRMANHTGKLDFFILMGLFRRSKHPALLSVVIFVVFLKALSGNAVLILLIHCD AHLH SPMYFFISQLSLMDMAYISVTVPKMLLDQVMGVNKKVSAPECGMQMFLYLTLAG SEFFLLA TMAYDRYVAICHPLRYPVLMNHRVCLFLASGCWFLGSDGFMLTPITMSFPFC RSWEIHH FFCEVPAVTILSCSDTSLYETLMYLCCVLMLLIPVTIISSSYLLILLTVHRMNSAEG RKK AFATCSSHLTVVILFYGAAVYTYMLPSSYHTPEKDMMVSVFYTLTPVLNPLIYS LRNKD VMGALKKMLTVRFVL
Source	in vitro E.coli expression system
Target Names	OR2T5
Protein Names	Recommended name: Olfactory receptor 2T5 Alternative name(s): Olfactory receptor OR1-62
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.