



Recombinant Human Olfactory receptor 56A3 (OR56A3)

Product Code	CSB-CF836771HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH54
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MTTHRNDTLSTEASDFLLNCFVRSPSWQHWLSLPLSLLFLLAVGANTTLLMTI WLEASLH QPLYLLSLLSLLDIVLCLTVIPKVLTIWFDFLRPISFPACFLQMYIMNCFLAMES CTFM VMAYDRYVAICHPLRYPsiITDHFVVKAAmFILTRNVLMTLPiPILSAQLRYCGRN VIEN CICANMSVSRlSCDDVTINhLYQFAGGWTLlGSDlLIFLSYTFILRAVLRlKAE GAVAK ALSTCGSHFmLILFFSTILLVfVLTHVAKKKVSPDVPVLLNVLHhVIPAAlnPIIYG VRT QEIKQGMQRLLKKGc</p>
Source	in vitro E.coli expression system
Target Names	OR56A3
Protein Names	Recommended name: Olfactory receptor 56A3 Alternative name(s): Olfactory receptor 56A6
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.