



Recombinant Human Olfactory receptor 2V2 (OR2V2)

Product Code	CSB-CF839403HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96R30
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METWVNQSYTDGFFLLGIFSHSTADLVLFSSVMAVFTVALCGNVLLIFLIYMDP HLHTPM YFFLSQLSLMDLMLVCTNVPKMAANFLSGRKSISFVGCGIQIGLFVCLVGSEGL LLGLMA YDRYVAISHPLHYPILMNQRVCLQITGSSWAFGIIDGLIQMVVVMNFPYCGLRK VNHFFC EMLSLLKLACVDTSLEKVFACCVFMLLPFSIIVASYAHILGTVLQMHSQAQAW KKALA TCSSHLTAVTLFYGAAMFIYLRPRHYRAPSHDKVASIFYTVLTPMLNPLIYSLRN REVMG ALRKGLDRCRIGSQH
Source	in vitro E.coli expression system
Target Names	OR2V2
Protein Names	Recommended name: Olfactory receptor 2V2 Alternative name(s): Olfactory receptor 2V3 Olfactory receptor OR5-3
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.