



Recombinant Human Olfactory receptor 10H5 (OR10H5)

Product Code	CSB-CF839883HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGA6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MQGLNHTSVSEFILVGFSAPHLQLMLFLLFLLMYLFTLLGNLLIMATVWSERSL HMPMY LFLCALSITEILYTVAIIPRMLADLLSTQRSIAFLACASQMFFSFSGFTHSFLLT V MGY DRYVAICHPLRYNVLMRLGCTCRVGCWAGGLVMGMVVTSAIFHLAFCGHK EIHHFFCH VPPLLKLACGDDVLVVAKGVLVCITALLGCFLILLSYAFIVAAILKIPSAEGRNK AFS TCASHLTVVVVHYGFASVIYLPKGPQSPEGDTLMGITYTVLTPFLSPIIFSLRN KELKV AMKKTCFTKLFPQNC
Source	in vitro E.coli expression system
Target Names	OR10H5
Protein Names	Recommended name: Olfactory receptor 10H5 Alternative name(s): Olfactory receptor OR19-25 Olfactory receptor OR19-26
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.