



Recombinant Human Olfactory receptor 8J3 (OR8J3)

Product Code	CSB-CF844060HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGG0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPENFTRVTEFILTGVSSCEPELQIPLFLVFLVLYVLTMAGNLGIITLTSVDSRLQ NPMY FFLRHLAIINLGNSTVIAPKMLMNFLVKKKTTSFYECATQLGGFLFFIVSEVMML AVMAY DRYVAICNPLLYMVVSRRLCLLLVSLTYLYGFSTAIVVSPCIFSVSYCSSNIINH FYCD IAPLLALSCSDTYIPETIVFISAATNLVFSMITVLVSYFNIVLSILRIRSPEGRKKAF ST CASHMIAVTVFYGTMLFMLYLPQQTNHSLDTPDKMASVFYTLVIPMLNPLIYSLRN NDVNVA LKKFMENPCYSFKSM
Source	in vitro E.coli expression system
Target Names	OR8J3
Protein Names	Recommended name: Olfactory receptor 8J3 Alternative name(s): Olfactory receptor OR11-173
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.