



Recombinant Human Olfactory receptor 10A4 (OR10A4)

Product Code	CSB-CF872436HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H209
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMWENWTIVSEFVLVSFSALSTELQALLFLLFLTIYLVTLMGNVLIILVTIADSAL QSPM YFFLRNLSFLEIGFNLVIVPKMLGLIIQD TTISFLGCATQMYFFFFFFGAAECCLL ATMA YDRYVAICDPLHYPVIMGHISCAQLAAASWFSGFSVATVQTTWIFSFPFCGPN RVNHFFC DSPPVIALVCADTSVFELEALTATVLFILFPFLLILGSYVRILSTIFRMPSAEGKHQ AFS TCSAHLVLSLFYSTAILTYFRPQSSASSESKLLSLSSTVVTPMLNPIIYSSRNK EVKA ALKRLIHRTLGSQKL
Source	in vitro E.coli expression system
Target Names	OR10A4
Protein Names	Recommended name: Olfactory receptor 10A4 Alternative name(s): HP2 Olfactory receptor-like protein JCG5
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.