



# Recombinant Human Olfactory receptor 10H2 (OR10H2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-CF016396HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O60403  |
| <b>Product Type</b>      | Transmembrane Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Sequence</b>          | MLGLNHTSMSEFILVGFSAPHLQLMLFLLFLLMYLFTLLGNLLIMATVWSERSL<br>HTPMY<br>LFLCVLSVSEILYTVAIIPRMLADLLSTQRSIAFLACASQMFFSFSFGFTHSFLLT<br>VMGY<br>DRYVAICHPLRYNVLMSPRGCACLVGCSWAGGSVMGMVV TSAIFQLTFCGSH<br>EIQHFLCH<br>VPPLLKLACGNNVPAVALGVGLVCIMALLGCFLILLSYAFIVADILKIPSAEGRN<br>KAFS<br>TCASHLIVVIVHYGFASVIYLPKPGPHSQEGDTLMATTYAVLTPFLSPIIFSLRNK<br>ELKV AMKRTFLSTLYSSGT |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>Target Names</b>      | OR10H2  |
| <b>Protein Names</b>     | Recommended name: Olfactory receptor 10H2 Alternative name(s): Olfactory receptor OR19-23   |
| <b>Expression Region</b> | 1-315   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |