



# Recombinant Mouse Taste receptor type 2 member 3 (Tas2r3)

<b>Product Code</b>	CSB-CF752648MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q7TQA7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	MFGFIEGVFLVLTITEFILGNLVNGFIVSINSSYWFKSKKISLSNFIITSLALFRIFLL W IIFIDSLIIVFSYQTHDSGIMMQLIDVFWTFTNHFSIWLISCLSVFYCLKIASFSHPS FL WLKWRASRVVVGMLWGALLLSCVSTMSLMNEFKIYSALTRSKDTPNMTEYIRL KRQEYNL MHVLGNLWKIPSLIVSLVAYLLLLLSLGKHTQQMQQYSIDSRDQSAEAHKRAM RIISSFL LFFLFYFLSFMILSSSRFLPETRIARIIGVVISMSYLVGDSFILVCNNKLKHTFVA MLP CECGHLKPGSKGPSAS
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Tas2r3
<b>Protein Names</b>	Recommended name: Taste receptor type 2 member 3 Short name= T2R3 Alternative name(s): T2R137 Short name= T2R37 mT2r41
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.