



Recombinant Human Putative olfactory receptor 2I1 (OR2I1P)

Product Code	CSB-CF847733HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGU4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLANQASAEERFLLLGFSDWPSLQPVLFAVLVLLCYLLTLTGNSALVLLAVRDPR LHTPMY YFLCHLALVDAGFTTSVVPPLLNLRGPALWLPRSHCTAQLCASLALGSAECV LLAVMAL DRAAAVCRPLRYAGLVSPRLCRTLASASWLSGLTNSVAQTALLAERPLCAPRL LDHFICE LPALLKLACGGDGDTTENQMFAARVVILLPPFAVILASYGAVARAVCCMRFGS GRRRAVG TCGSHLTAVCLFYGSAIYTYLQPAQRYNQARGKFVSLFYTVVTPALNPLIYTLR NKKVKG AARRLLRSLGRGQAGQ
Source	in vitro E.coli expression system
Target Names	OR2I1P
Protein Names	Recommended name: Putative olfactory receptor 2I1 Alternative name(s): Putative olfactory receptor 2I2 Putative olfactory receptor 2I3 Putative olfactory receptor 2I4
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.