



Recombinant Human Olfactory receptor 8J1 (OR8J1)

Product Code	CSB-CF850897HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGP2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPENFTRVTEFILTGVSSCPQLIPLFLVFLVLYGLTMAGNLGIITLTSVDSRLQ TPMY FFLQHLALINLGNSTVIAPKMLINFLVKKKTTSFYECATQLGGFLFFIVSEVIMLAL MAY DRYVAICNPLLYMVVSRRLCLLLVSLTYLYGFSTAIVVSSYVFSVSYCSSNIINH FYCD NVPLLALSCSDTYLPETVVFISAATNVVGSLLIIVLVSYFNIVLSILKICSSEGRKKA FST CASHMMAVTIFYGTLLFMVYQPRSNHSLDTHDDKMASVFYTLVIPMLNPLIYSLR NKDVKT ALQRFMTNLCYSFKTM
Source	in vitro E.coli expression system
Target Names	OR8J1
Protein Names	Recommended name: Olfactory receptor 8J1 Alternative name(s): Olfactory receptor OR11-183
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.