



# Recombinant Human Olfactory receptor 10H4 (OR10H4)

|                          |                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF854125HU                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                               |
| <b>Uniprot No.</b>       | Q8NGA5                                                                                                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Transmembrane Protein                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                            |
| <b>Sequence</b>          | MPSQNYSIIEFNLFGFSAFPQHLLPILFLLYLLMFLFTLLGNLLIMATIWIEHRLH<br>TPM<br>YLFLCTLSVSEILFTVAITPRMLADLLSTHHSITFVACANQMFFSFMFGFTHSFL<br>LVMG<br>YDRYVAICHPLRYNVLMSPRDCAHLVACTWAGGSVMGMMVTTIVFHLTFCGS<br>NVIHHFFC<br>HVLSELLKLACENKTSSVIMGVMLVCVTALIGCLFLIILSYVFIVAAILRIPSAEGRH<br>KTF<br>STCVSHLTVVVTHYSFASFIYLPKGLHSMYSDALMATTYTVFTPFLSPIIFSLR<br>NKELK NAINKNFYRKFCPPSS |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | OR10H4                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: Olfactory receptor 10H4 Alternative name(s): Olfactory receptor OR19-28                                                                                                                                                                                                                                                                       |
| <b>Expression Region</b> | 1-316                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                   |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                              |