



Recombinant Human Olfactory receptor 12D3 (OR12D3)

Product Code	CSB-CF887020HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UGF7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MENVTTMNEFLLLGLTGVQELQPFFFGIFLIYLINLIGNGSILVMVVLEPQLHSP MYFF LGNLSCLDISYSSVTLPKLLVNLVCSRRAISFLGCITQLHFFHFLGSTEAILLAIMA FDR FVAICNPLRYTVIMNPQVCILLAAAALWISFFYALMHSVMTAHLSCGSGQKLNH FFYDVK PLLELACSDTLNQQWLLSIVTGSISMGAFFLTLLSCFYVIGFLLFKNRSCRILHKA LSTC ASHFMVVCLFYGPVGFYIRPASATSMIQDRIMAIMYSAVTPVLNPLIYTLRNKE VMMAL KKIFGRKLFKDWQQHH
Source	in vitro E.coli expression system
Target Names	OR12D3
Protein Names	Recommended name: Olfactory receptor 12D3 Alternative name(s): Hs6M1-27 Olfactory receptor OR6-27
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.