



Recombinant Human Taste receptor type 2 member 3 (TAS2R3)

Product Code	CSB-CF889154HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NYW6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMGLTEGVFLILSGTQFTLGILVNCFIELVNGSSWFKTKRMSLSDFIITTLALLRII LLC IILTDSFLIEFSPNTHDSGIIMQIIDVSWTFTNHLSIWLATCLGVLYCLKIASFSHPT FL WLKWRVSRVMVWMLL GALLSCGSTASLINEFKLYSVFRGIEATRNVTEHFRK KRSEYYL IHVLGTLWYLPPLIVSLASYSLLIFSLGRHTRQMLQNGTSSRDPTTEAHKRAIRII LSFF FLFLLYFLAFLIASFGNFLPKTKMAKMIGEVMTMFYFAGHSFILILGNSKCLKQTFV VMLR CESGHLKPGSKGPIFS
Source	in vitro E.coli expression system
Target Names	TAS2R3
Protein Names	Recommended name: Taste receptor type 2 member 3 Short name= T2R3
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.