



Recombinant Pan troglodytes Taste receptor type 2 member 14 (TAS2R14)

Product Code	CSB-CF737380EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q646B7
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MGGVIKSIFTFVLIVEFIIGNLGN SFIALVNCIDWVKGRKISSVDRILTALAIKISLV W LIFGSWCVSVFFPALFATEKMFRLTNIWTVINHFSVWLATGLGTFYFLKIANFS NSIFL YLKWRVKKVVLVLLLVT SVFLFLNIALINIHINASINGYRRNKTCSSDSSNFTRFS SLIV LTSTVFIFIPFTLSLAMFLLLIFSXWKHRKKMQHTVKRSGDASTKAHRGVKSVXT FFLLY AIFCLSFFISVWTSERLEENLIILSQVMGMAYPSCHSCVLILGNKKLRQASLSVLL WLRV MFKDGEP SGHKEFRESS
Source	in vitro E.coli expression system
Target Names	TAS2R14
Protein Names	Recommended name: Taste receptor type 2 member 14 Short name= T2R14
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.