



Recombinant Mouse Olfactory receptor Olfr180 (Olfr180)

Product Code	CSB-CF745881MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TS48
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEKTNHSLTTQFILVGFSDHPDLKTPLFLLFSVIYLVTMVGNLGLVAVIYLEPRL HTPMY IFLGNLALMDSCCSCAITPKILENFFSVDRRISLYECMAQFYFLCLAETADCFL AAMAY DRYVAICNPLQYHSMMSKKLSIQMSIGTFITSNLHSLIHVGCLLRLTFCKSNRID HFFCD ILPLYRLSCTDPFINELMIYIFSMPIQVFTITTVLVSYFCILLTIFKMKSKDGRGKAF ST CASHFFSVSIFYVCLLMYIRPFDEGNKDIPVAVFYTIIPLLNPFYISLRNKEVVNA VKK VMKTHSIFKNASASMAR
Source	in vitro E.coli expression system
Target Names	Olfr180
Protein Names	Recommended name: Olfactory receptor Olfr180
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.