



Recombinant *Ashbya gossypii* ADIPOR-like receptor IZH1 (IZH1)

Product Code	CSB-CF748344DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q753H5
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	MSEERGMKEQTISEMAQEMAHTTSEGLKKRIRKLYTFDELPAWQKDNELILSG YVRETNS VKECLRAMTYFNNE SINIYTHLIPGVAYLVFLIFADLVLAQLLPGLDAGEHRML RFYLL GAFTCLACSSCFHCLKQHSEPHSRLWSKVDYLGILAQITCSTISLLYYGYHSYP SHFVFF STLTVALCSACAVLVLND SFNTVAFRPLRAFLFMAFGLSGVIPVLAGSYQFGFA EWAARI QLKYVLYEAVFYITGALVYGFRIPERFAPGKFD MVGHSHQIFHLLVVLGTLCHF RAVTGS YIFICTGKHYSSLLMFI
Source	in vitro E.coli expression system
Target Names	IZH1
Protein Names	Recommended name: ADIPOR-like receptor IZH1
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.