



# Recombinant Cercopithecus neglectus Melanocyte-stimulating hormone receptor (MC1R)

<b>Product Code</b>	CSB-CF769680CGW
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q864K1
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Cercopithecus neglectus (De Brazza's monkey)
<b>Sequence</b>	MPVQGSQRLLGSLNSTPTAPPHLGLAANQTGTRCLEVSIPDGLFSLGLVSL VENLVV TAIAKNRNLHSPMYCFICCLALSDLLVSGSNMLETAVILLEAGALAARAADVQQ LDNVI DVITCSSMLSSLCLFLGAIIVDRYISIFYALRYHSIVTLPRARRAVAAIWWASVLFS MLFI AYYDHA AVL LCLV VFFLAMLV LMAVLYIHMLARACQHAQGIARLHKRQCPAHQ GFGLKGA ATLTILLGIFFLCWGPFFLH LTLIVLCPQHPTCSCIFKNFNLFLALIICNAIIDPLIYA F RSQELRRTLKEVLLCSW
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MC1R
<b>Protein Names</b>	Recommended name: Melanocyte-stimulating hormone receptor Short name= MSH-R Alternative name(s): Melanocortin receptor 1 Short name= MC1-R
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.