



Recombinant Human Olfactory receptor 10AD1 (OR10AD1)

Product Code	CSB-CF844058HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGE0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLRNGSIVTEFILVGFQQSSTSTRALLFALFLALYSLTMAMNGLIIFITSWTDPKL NSPM YFFLGHLSSLVDVCFITTTIPQMLIHLVVRDHIVSFVCCMTQMYFVFCVGVAECILL AFMA YDRYVAICYPLNYVPIISQKVCVRLVGTAWFFGLINGIFLEYISFREPFRRDNHIE SFFC EAPIVIGLSCGDPQFSLWAIFADAIVVILSPMVLTVTSYVHILATILSKASSSGRG KTFS TCASHLTVVIFLYTSAMFSYMNPHSTHGPDKDKPFSLLYTIITPMCNPIIYSFRN KEIKE AMVRALGRTRLAQPQSV
Source	in vitro E.coli expression system
Target Names	OR10AD1
Protein Names	Recommended name: Olfactory receptor 10AD1 Alternative name(s): Olfactory receptor OR12-1
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.