



Recombinant Gorilla gorilla gorilla Taste receptor type 2 member 60 (TAS2R60)

Product Code	CSB-CF731062GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q645Z0
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	MNGDHMVLGSSVTDKKAAILVTILLLLRLVAIAGNGFIIAALGVEWVLRRLMLPCD XLLV SLGASRFCLQSVVMGKTIYVFLHPMAFPYNPVLQFLAFQWDFLNAATLWFST WLSVFYCV KIAAFTHPVFLWLKHKLSGWLPWILFSSVGLSSFTTILFFIGNHRMYQNYLRNHL QPWNI TGNSIRSYCEKFYLFPLKMITWTMPTAVFFICMILLITSLGRHMKKALLTTSGFR EPSMQ AHIKALLALLSFAMLFISYFLSLVFSAAGIFPPLDFKFWVWVESVIYLCAAVHPIILL FSN CRLRAVLKSCRSSRCGTP
Source	in vitro E.coli expression system
Target Names	TAS2R60
Protein Names	Recommended name: Taste receptor type 2 member 60 Short name= T2R60 Alternative name(s): T2R56
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.