



Recombinant Human Olfactory receptor 1B1 (OR1B1)

Product Code	CSB-CF850899HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGR6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMSFAPNASHSPVFLLLGFSRANISYTLFFFLFLAIYLTTLGNVTLVLLISWDSR LHSP MYLLRGLSVIDMGLSTVTLPLQLLAHLVSHYPTIPAARCLAQFFFFYAFGVTDTL VIAVM ALDRYVAICDPLHYALVMNHQRCACLLALSWWVSILHTMLRVGLVPLCWTGD AGGNVNL PHFFCDHRPLL RASCSDIHSNELAIFFEGLMLGPCALIVLSYVRIGAAAILRLPS AAGR RRAVSTCGSHLTMVGFLYGTIICVYFQPPFQNSQYQDMVASVMYTAITPLANP FVYSLHN KDVKGALCRLLEWVKVDP
Source	in vitro E.coli expression system
Target Names	OR1B1
Protein Names	Recommended name: Olfactory receptor 1B1 Alternative name(s): Olfactory receptor 9-B Short name= OR9-B Olfactory receptor OR9-26
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.