



Recombinant *Oenothera biennis* Cytochrome c biogenesis protein ccsA (ccsA)

Product Code	CSB-CF539065ODZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0Z508
Product Type	Transmembrane Protein
Immunogen Species	<i>Oenothera biennis</i> (German evening primrose) (<i>Onagra biennis</i>)
Sequence	MIFYTLEHILTHISFSLVSIGITIFLITLSVDEIIGLYDSSEKGVIGTFLCITGLLVTRW AYSGHFPLSNLYESLLFLSWSFAIIHMFYPYFKKQKSYVRTITSSSTIFTQGLVTS GLLSE MQQSEILVPALQSQWLMMHVSMMVLGYAALLCGSLLSVALLVITFRKALRIFSK KKAFLK DSFSFVEIQYRNEPSNVLLSTSFISSKNYYRAQLIQQLDRWSSRIISLGFILTIGI LSG AVWANEAWGSYWNWDPKETWAFITWTMFAIYLHTRTNPNFQSVNSAIVAFLG FIIIWICY FGVNLLGIGLHSYGSFNLH
Source	in vitro E.coli expression system
Target Names	ccsA
Protein Names	Recommended name: Cytochrome c biogenesis protein ccsA
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.