



Recombinant Bovine Probable G-protein coupled receptor 171 (GPR171)

Product Code	CSB-CF673604BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3ZBK9
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MTNSSTFCPVYRDLEPFTYFFYLVLFIGIIGSCFATWAFIQKNTNHRCVSIYLINL LTAD FLLTLALPVKITVDLGVAPWKLRIHFHCQVTACLIYINMYLSIIFLAFVSDRCLQLTY SC KIYRIQEPGFAMISAVVWLMVLLIMVPNMIPIKDIKEKPNVGCMEFKSEFGRN WHLLT NFISIAIFFNFSAIILISNCLVIRQLYRNKDNENYPNVKRALISILLVTTGYIICFVPY H IVRIPYTLSQTEVISDCSTRISLKFKAKEATLLAVSNLCFDPILYYHLSKAFRLKITE TF ASHKESKAQKEKPRSENNNA
Source	in vitro E.coli expression system
Target Names	GPR171
Protein Names	Recommended name: Probable G-protein coupled receptor 171
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.