



Recombinant Pongo pygmaeus Taste receptor type 2 member 30 (TAS2R30)

Product Code	CSB-CF731061EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q645V4
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	MITFLPIIFSILIVVIFVIGNFANGFIALVNSIEWVKRQKISFADQILIALAVSRVGLL W ALLLHWYATELNLAFYSEVRITAYNVWAVTNHFSNWLATSLSMFYLLKIANFS NLIFLR IKRRVKSVILVILLGPLLFLVCHLFVINMNEIVWTKEYEGNLTWKIKLRNAVFLSN MTLT MLANFVPLTLTLISFLLLICSLCKHLKMKQLHGKGSQDPSTKVHIKALQTVTCFL LLCAI YFLSMIISVYNFGRLEKKPVFMFCQAITFSYPSTHAFILIWGNKKLKQIFLSVLWH VRYW VKDRSLRLHRFTRAALCKG
Source	in vitro E.coli expression system
Target Names	TAS2R30
Protein Names	Recommended name: Taste receptor type 2 member 30 Alternative name(s): Taste receptor type 2 member 47 Short name= T2R47
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.