



Recombinant Mouse Taste receptor type 2 member 39 (Tas2r39)

Product Code	CSB-CF766512MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TQA5
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAQPSNYWKQDVLPLSILMLTLVATECTIGIIASGIVMAVNAVSWVQKKAISITT RILLL LSVSRIGLQSIMLIEITSSIFNVAFYNSVLYRVSNSVSVFLNYCSLWFAALLSFFH FVKI ANFSYPLFFKLKWRISELMPWLLWLSVFISFSSSMFFSKHKFTVNNNNSLSNNI CNFTMK LYVVETNVVNSFLFISGILPPLTMFVATATLLIFSLRRHTLNMRNSATGSRNPCI EAHM QAIKETSCFLFLYILNAAALLLSTSNIVDASLFWISIVIRIVLPVYPAGHSVLLIQNN PGL RRTWKHLQSQIHLYLQNR
Source	in vitro E.coli expression system
Target Names	Tas2r39
Protein Names	Recommended name: Taste receptor type 2 member 39 Short name= T2R39 Short name= mT2R34
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.