



Recombinant Rabbit Adenosine receptor A3 (ADORA3)

Product Code	CSB-CF001379RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02667
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MPDNSTTLFLAIRASYIVFEIVIGVCAVVGNVLVIWVIKLNPSLKTTFYFIFSLAL ADI AVGFLVMPLAIVISLGITIGFYSCLVMSCLLLVFTHASIMSLLAIVDRYLRVKLTV RYR RVTTQRRIWALGLCWVVSLLVGFTPMFGWNMKPTLESARNYSDFQCKFDSV IPMEYMFV FSFFTWILIPLLLMCALYVYIFYIIRNKLVSFSSFKETGAFYRREFKTAKSLFLVL ALF AGCWLPLSIINCVTYFKCKVPDVVLLVGILLSHANSMMNPVYACKIQKFKETYL LIFKA RVTCQPSDSLDPSSSEQNSE
Source	in vitro E.coli expression system
Target Names	ADORA3
Protein Names	Recommended name: Adenosine receptor A3
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.