



Recombinant Human Olfactory receptor 2T12 (OR2T12)

Product Code	CSB-CF850883HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NG77
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEMRNTTPDFILLGLFNHTRAHQVLFMMLLATVLTSLFSNALMILLIHWDHRLH RPMYFL LSQLSLMDMMLVSTTVPKMAADYLTGNKAISRAGCGVQIFFLPTLGGGECFLL AAMAYDR YAAVCHPLRYPTLMSWQLCLRMTMSSWLLGAADGLLQAVATLSFPYCGAHEI DHFFCEAP VLVRLACADTSVFENAMYICCVLMLLVPFSLILSSYGLILAAVLLMRSTEARKKA FATCS SHVAVVGLFYGAGIFTYMRPKSHRSTNHDKVVSIFYTMFTPLLNPLIYSVRNSE VKEALK RWLGTCVNLKHQQNEAHRSR
Source	in vitro E.coli expression system
Target Names	OR2T12
Protein Names	Recommended name: Olfactory receptor 2T12 Alternative name(s): Olfactory receptor OR1-57
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.