



Recombinant Human Olfactory receptor 7C1 (OR7C1)

Product Code	CSB-CF017055HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O76099
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METGNQTHAQEFLLLGFSATSEIQFILFGLFLSMYLVTFGTGNLLIILAICSDSHLH TPMY FFLSNLSFADLCFTSTTVPKMLLNILTQNKFITAGCLSQIFFFTSFGCLDNLLLT VMAY DRFVAVCHPLHYTVIMNPQLCGLLVLGSWCISVMGSLLETTLVLRLSFCTEME I PHFFCD LLEVLKLACSDTFINNVIYFATGVLGVISFTGIFFSYKIVFSILRISSAGRKHKAF ST CGSHLSVVTLFYGTGFGVYLSSAATPSSRTSLVASVMYTMVTPMLNPFIYSLR NTDMKRA LGRLLSRATFFNGDITAGLS
Source	in vitro E.coli expression system
Target Names	OR7C1
Protein Names	Recommended name: Olfactory receptor 7C1 Alternative name(s): Olfactory receptor 7C4 Olfactory receptor OR19-16 Olfactory receptor TPCR86
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.