



Recombinant Rat Adenosine receptor A3 (Adora3)

Product Code	CSB-CF001379RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28647
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MKANNTTTSALWLQITYITMEAAIGLCAVVGNNMLVIWVVKLNRTLRTTTFYFIVS LALAD IAVGVLVIPLAIAVSLEVQMHPYACLFMSCVLLVFTHASIMSLAIAVDRLRVKL TVRY RTVTTQRRIWFLGLCWLVSFVGLTPMFGWNRKVTLELSQNSSTLSCHFRRSV VGLDYMV FFSFITWILIPLVVMCIYLDIFYIIRNKLSQNLTFRETRAFYGREFKTAKSLFLVL FL FALCWLPISIINFVSYFNVKIPEIAMCLGILLSHANSMMNPVYACKIKKFKETYF VILR ACRLCQTSDSLDSNLEQTTE
Source	in vitro E.coli expression system
Target Names	Adora3
Protein Names	Recommended name: Adenosine receptor A3 Alternative name(s): TGPCR1
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.