




Recombinant Human Olfactory receptor 5V1 (OR5V1)

Product Code	CSB-CF871401HU
Abbreviation	Recombinant Human OR5V1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UGF6
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	MERKNQTAITEFIILGFSNLNELQFLLFTIFFLTYFCTLGGNIIILTTVTDPHLHTP MY YFLGNLAFIDICYTTSNVPQMMVHLLSKKKSISYVGCVVQLFAFVFFVGSECLLL AAMAY DRYIAICNPLRYSVILSKVLCNQLAASCWAAGFLNSVHTVLTFCLPFCGNNQI NYFFCD IPLLILSCGNTSVNELALLSTGVFIGWTPFLCIVLSYICIIISTILRIQSSEGRRKAF ST CASHLAIVFLFYGSAIFTYVRPISTYSLKKDRLVSVLYSVVTPMLNPIIYTLRNKDI KEA VKTIGSKWQPPISSLD SKLTY
Source	in vitro E.coli expression system
Target Names	OR5V1
Protein Names	Recommended name: Olfactory receptor 5V1 Alternative name(s): Hs6M1-21 Olfactory receptor OR6-26
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.