



Recombinant Human Olfactory receptor 3A3 (OR3A3)

Product Code	CSB-CF016638HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47888
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSLQKLMEPEAGTNRTAVAEFILLGLVQTEEMQPVVFVLLLFAYLVTGGNLSI LAAVLV EPKLHAPMYFFLGNLVLDVGCITVTPAMLGRLLSHKSTISYDACLSQLFFFH LLAGMD CFLLTAMAYDRLLAICQPLTYSTRMSQTVQRMVLAASWACAFTNALHTVAMS TLNFCGP NEVNHFYCDLPQLFQLSCSSTQLNELLLFVAAAFMAVAPLVFISVSYAHVVAAV LQIRSA EGRKKAFASTCGSHLTVVGIFYGTGVFSYMLGSGVESDSDKDGVGVFMTVINP MLNPLIYS LRNTDVQGALCQLLVGKRSLT
Source	in vitro E.coli expression system
Target Names	OR3A3
Protein Names	Recommended name: Olfactory receptor 3A3 Alternative name(s): Olfactory receptor 17-201 Short name= OR17-201 Olfactory receptor 3A6 Olfactory receptor 3A7 Olfactory receptor 3A8 Olfactory receptor OR17-22
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.