



Recombinant Human Olfactory receptor 5K4 (OR5K4)

Product Code	CSB-CF016964HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NMS3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MARENHSLAAEFILIGFTNYPELKTLFFVVFSAIYLVTMVGNLGLVALIYVERRLL TPMY IFLGNLALMDSCCSCAVTPKMLENFFSEDRIISLYECMAQFYFLCLAETTDCFL ATMAY DRYVAICHPLQYHTMMSKTL CIRMTTGAFKAGNLHSMIHVGLLLRLTFCRSNKI HHFFCD ILPLYRLSCTDPSINELMIYIFSIPIQIFTIATVLISYLCILLTVFKMKSKEGRGKA FS T CASHFLSVSIFYICLLMYIGPSEEGDKDTPVAIFYAIVIPLLNPFIYSLRNKEVIN VLL KK IMRNYNILKQTCSIANLFLIY
Source	in vitro E.coli expression system
Target Names	OR5K4
Protein Names	Recommended name: Olfactory receptor 5K4
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.