



Recombinant Mouse Mas-related G-protein coupled receptor member X1 (Mrgprx1)

Product Code	CSB-CF817224MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8CIP3
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MDPTISSHDTESTPLNETGHPNCTPILTLSFLVLITTLVGLAGNTIVLWLLGFRM RRKAI SVYILNLALADSFLLCCHFIDSLRLRIIDFYGLYAHKLSKDILGNAAIIPYISGLSILSA I STERCLCVLWPIWYHCHRPRNMSAICALIWVLSFLMGILDWFSGFLGETHHHL WKNVDF IITAFILFLMLLSGSSLALLLRILCGPRRKPLSRLYVTIALTVMVYLICGLPLGLYL FL LYWFGVHLHYPFCHIQVTAVLSCVNSSANPIIYFLVGSFRQHRKHRSLKRVLK RALEDT PEEDEYTDSHLHKTTEISESRY
Source	in vitro E.coli expression system
Target Names	Mrgprx1
Protein Names	Recommended name: Mas-related G-protein coupled receptor member X1 Alternative name(s): Mas-related G-protein coupled receptor member C11 Sensory neuron-specific G-protein coupled receptor 1
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.