



Recombinant Human Olfactory receptor 1E2 (OR1E2)

Product Code	CSB-CF016495HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47887
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMGQNQTSISDFLLLGLPIQPEQQNLCYALFLAMYLTLLGNLLIIVLIRLDShLH TPVY LFLSNLSFSDLCFSSVTMPKLLQNMQNQDPSIPYADCLTQMYFFLYFSDLESFL LVAMAY DRYVAICFPMHYTAICFLLHYTAIMSPMLCLSVVALSWVLTTFHAMLHTLLMAR LCFCAD NVIPHFFCDMSALLKLACSDTRVNEWVIFIMGGLILVIPFLLILGSYARIVSSILKV PSS KGICKAFSTCGSHLSVVSIFYGTVIGLYLCPSANSSTLKDTVMAMMYTVVTPML TPFIYS LRNRDMKGALERVICKRKNPFL
Source	in vitro E.coli expression system
Target Names	OR1E2
Protein Names	Recommended name: Olfactory receptor 1E2 Alternative name(s): Olfactory receptor 17-93/17-135/17-136 Short name= OR17-135 Short name= OR17-136 Short name= OR17-93 Olfactory receptor 1E4
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.