



Recombinant Human Olfactory receptor 1L3 (OR1L3)

Product Code	CSB-CF844082HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH93
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGMSNLTRLSEFILLGLSSRSEDQRPLFALFLIYLVTLMGNLLILAIHSDPRLQN PMY FFLSILSFADICYTTVIVPKMLVNFLSEKKTISYAECLAQMYFFLVFGNIDSYLLAA MAI NRCVAICNPFHYVTVMNRRCCVLLAFPITFSYFHSLLHVLLVNRLTFCTSNVIH HFFCD VNPVLKLSCSSTFVNEIVAMTEGLASVMAPFVCIISYLRILIAVLKIPSAAGKHKA FST CSSHLTVILFYGSISYVYLQPLSSYTVKDRIATINYTVLTSVLNPFYSLRNKDM KRGL QKLINKIKSQMSRFSTKTNKICGP
Source	in vitro E.coli expression system
Target Names	OR1L3
Protein Names	Recommended name: Olfactory receptor 1L3 Alternative name(s): Olfactory receptor 9-D Short name= OR9-D Olfactory receptor OR9-28
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.